

HD TrueChrome IIS



	TrueChrome II	TureChrome IIS	HD Lite
Advantage	Perfect color reproduction Highly intelligent	Calibration and Measurement Perfect color reproduction	5MP
	High speed data transmission Fluorescence imaging		
Sensor	IMX236		MT9P006
Sensor size	1/2.7"		1/2.5"
Dynamic resolution	1920x1080		
Static resolution	3264x1836		
Frame rate	1920x1080 30fps via USB2.0 1920x1080 60fps via HDMI		1920x1080 15fps via USB2.0 1920x1080 15fps via HDMI
Capture	High speed SD card(8G)		High speed SD card(8G)
Video recording	1080p 30fps @SD Card 1080p 30fps @PC		1080p 15fps @SD Card 1080p 15fps @PC
Exposure time	0.001sec-10.0sec		
White balance	Auto		
Setting	Gain , gamma , saturation , contrast		
Built-in software	Cloud 1.0 Ver		
PC software	ISCapture		
Output port	USB2.0、microSD Card、HDMI		
Compatible OS	Window XP and above , MAC OSX		
Optical port	Standard C-Mount		

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt P.O.Box: 127 Giza, Egypt Code:12515

Tel. Fax:+(202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)

www.photoscience-co.com / info@photoscience-co.com