

AUTOMATIC WATER STILL

Use : Automatic Water Still produces and reserves distilled water automatically.

This unit offers exceptional value for the user. It is very easy and safe to operate with indicators such as reservoir full indicator, distilling indicator, and with safety device such as overheat safety sensor



Inside of the unit




Water still dispenser



Auto sensor for tank full



Filter



Valve




Overheat prevention device



Specification

Model NO	PHAWS-FA12		PHAWS-FA30	
Distilled Water Capacity	3.5~4.0l/ hour		3.5~4.0l/ hour	
Dimensions (WxDxH,mm)	Boiling Water Tank	200x200x300	200x200x300	
	Distilled Water Tank	200x200x300	275x275x400	
	Out size	530x300x590	680x375x600	
Heater	3.0 kW			
Reservoir capacity	Approx. 8l		Approx. 20l	
Material	Internal	Stainless Steel(#304)		
	External	Steel, Powder coating		
Safety Device	Over temp limiter.			
	Water Supply Cut-Off			
Operation	Fully Automatic Operation			
Signal Lamp	Distillation, Filled with water Lamp			
Electric Supply	220V, 50/60 Hz			



Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt P.O.Box: 127 Giza, Egypt Code:12515
Tel. Fax: +(202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)
www.photoscience-co.com / info@photoscience-co.com