

# Balance FA-B(0.1mg/0.0001g)



FA-B: Electromagnetic  
 Balance with RS232 interface.  
 External calibration, easy calibration operating.  
 Big glass wind shield make sure the reading stable.  
 We also provide OEM service,

## Specification Series FA-B

Model	FA1004B	FA1104B	FA1604B	FA2004B	FA2204B
Max Cap.	100g	110g	160g	200g	220g
Readability	0.1mg				
Pan Size(mm)	Φ80				
Calibration	External Calibration				
Display	LCD(White Back Light)				
I/O	RS232				
Power	220V/110V				
Unit	g, ct, oz, pcs				
Standard	RS232, wind shield, operating manual, Calibration weight				
Optional	Printer(dot-matrix printer, thermal printer)				
Outline Dimension(cm)	34.5×21.5×34				
Packing Size(cm)	39×30×45				
G.W.(kg)	10				

**Standard**
- ✓ +  
TOL
🎯  
LEVEL
0.00  
TARE
1 YEARS  
GUARANTY
🧪  
CAL EXT
🛡️  
SHIELD
📡  
RS 232
🏠  
PCS
%  
PERCENT
🏭  
100V-240V
♻️  
UNIT

**Optional**
0.00G  
0.00G  
Dual Display
🖨️  
UNDER

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515  
**Tel. Fax:** +(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)  
[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)

# Balance FA-A(0.1mg/0.0001g)



**FA-A series:**

**Electromagnetic balance. External calibration,  
Easy calibration operating. RS232 interface,  
Big glass wind shield make sure the reading is stable.  
We also provide OEM service,**

## Specification Series FA-A

Model	FA1004A	FA1104A	FA1604A	FA2004A	FA2204A
Max Cap.	100g	110g	160g	200g	220g
Readability	0.1mg				
Pan Size(mm)	Φ80				
Display	LCD(White Back Light)				
I/O	RS232				
Power	220V/110V				
Unit	g, ct, oz, ozt, %, pcs				
Standard	Calibration <u>weight</u> ,RS232				
Optional	printer(dot-matrix printer, thermal printer)				
Outline Dimension(cm)	35×21×35				
Packing Size(cm)	39×30×45				
G.W.(kg)	8.5				

**Standard**

TOL	LEVEL	TARE	1 YEARS GUARANTY	CAL EXT	SHIELD	RS 232	PCS	PERCENT	100V-240V	UNIT
-----	-------	------	------------------	---------	--------	--------	-----	---------	-----------	------

**Optional**


MOVE	Dual Display	UNDER
------	--------------	-------

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515  
**Tel. Fax: + (202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)**

[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)

# Balance FA-C(0.1mg/0.0001g)

FA-C series,  
 Electromagnetic balance. 0.1mg accuracy,  
 Max range: 220g. Internal calibration,  
 With RS232 interface.  
 Big glass wind shield make sure the reading is stable.  
 We also provide OEM service,



## Specification Series FA-C

Model	FA1004C	FA1104C	FA1604C	FA2004C	FA2204C
Max Cap.	100g	110g	160g	200g	220g
Readability	0.1mg				
Pan Size(mm)	Φ80			Φ110	
Calibration	Automatically				
Display	LCD(White Back Light)				
I/O	RS232				
Power	220V/110V				
Unit	g, ct, oz, pcs				
Standard	RS232, Internal calibration, wind shield, operating manual				
Optional	printer ( <u>dot-matrix printer</u> , <u>thermal printer</u> )				
Outline Dimension(cm)	34.5×21.5×34				
Packing Size(cm)	39×30×45				
G.W.(kg)	10				

**Standard**

TOL	LEVEL	TARE	1 YEARS GUARANTY	CAL EXT	CAL INT	SHIELD	RS 232	PCS	PERCENT	100V-240V	UNIT


**Optional**

Dual Display	UNDER

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515  
**Tel. Fax:** +(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)  
[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)

# Balance FA-J(0.1mg/0.0001g)

FA-J series,  
 Special gravity density balance for solid and liquid.  
 0.1mg accuracy,  
 With RS232 interface.  
 Easy operating,  
 We also provide you OEM service,



## Specification Series FA-J

Model:	FA2104J	FA2004J	FA1004J	JA5003J	JA3003J	JA2003J
Max. Cap.	210g	200g	100g	500g	300g	200g
Readability:	0.1mg			1mg		
Pan Size(mm):	Φ80mm~90mm					
Display	LCD					
I/O:	RS232					
Power:	220V/110V					
Unit	ct,oz,etc					
Standard	Density kits;calibration <u>weight</u>					
Optional						
Outline Dimension(cm)	34.5×21.5×34					
Packing Size(cm):	48×38×52					
G.W.(kg):	10kg					

Standard

--	--	--	--	--	--	--	--	--	--	--

Optional

--	--

Dual Display

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515  
**Tel. Fax:** +(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)  
[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)