

## BL series portable Gem microscope

EXW USD 35.-/pc

Small size mini pocket Gem microscope

New design mini powerful darkfield stage

The best travel Gem testing equipment

### Main specifications

Eyepieces: WF10x or optional eyepieces

Focus system: Eyepieces with rotation focusing

Stage: Mini darkfield stage can be added Gem stone clips beside

Illumination: LED illumination



## LED illumination mini pocket polariscope

EXW USD 35.-/pc

May be the smallest mini pocket polariscope

Easy to add a Gem stone clips beside

With LED illumination

### Main specifications

Polarizers: Rotation knob on upper one polarized lens

Stage: Simple stage with Gem stone clips beside

Illumination: LED illumination



## BL series portable microscope

EXW USD 28.-/pc(WF10x eyepieces, 10x objective, natural lighting)

EXW USD 35.-/pc(WF10x eyepieces, 10x objective, transmitted illumination)

EXW USD 48.-/pc(WF10x eyepieces, 4x,10x,40x objective, transmitted illumination)

Optional projection illumination unit EXW USD 5.-/pc

Small size instant use Mini pocket handheld portable microscope, multipurpose microscope

Can work as biological microscope or simple stereo microscope

## Main specifications

Eyepieces: WF10x DIN18 or optional eyepieces

Objective: 10x, optional 4x and 40x

Focus system: Rotation focusing knob, coarse focusing ring on main body ring on body of microscope, fine focusing on slide stage

Stage: inserted stage to accept a slide

Illumination: Transmitted LED illumination base, optional projection illumination

## View with natural lighting

## Slide stage and transmitted LED illumination base



## Mini pocket handheld Gem microscope, natural lighting

Small size mini pocket Gem microscope

The best travel Gem testing equipment

## Main specifications

Eyepieces: WF10x DIN18 or optional eyepieces

Objective: 1x or customized magnification

Focus system: Rotation focusing knob

Gem stone clips: New design Gem stone clips, simple use

Illumination: Natural light

