

INCUBATOR SHAKER

Use This product is used for cultivation and shaking test of circle movement .



**Safety device
(Stop for opening)**




Magnet packing



The shelves easy to move



PID controller



Spring wire rack



Features

- Accurate temperature controls by Digital PID controller.
- Accurate speed control and digital mark function made for user convenience.
- Circulation device in internal chamber have a little noise and keep overall temperature evenly.
- Overheat, over current prevention device for users and safety function to stop shaking device automatically when door open.
- Various holder and plate available according to users spec.
- Little noise designed structurally.

Model NO		PH-SI64L	PH-SI150L
Temp	Range	Amb +5 °C ~ 60 °C	
	Controller	Digital PID Controller	
	Accuracy	±0.2 °C	
Size	Internal(WxDxH, mm)	400x400x400	560x560x500
	External(WxDxH, mm)	470x580x630	630x740x750
Safety	Electrical	Over temp limiter, Door safety	
Material	Internal	Stainless steel(STS #304)	
	External	Steel, Powder coating	
Speed Range		20 ~300 rpm	
Shaking Stroke		20 mm Orbital Motion	
Electric Requirements		220VAC. 50/60Hz	

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt / **P.O.Box:** 127 Giza, Egypt **Code:** 12515

Tel. Fax: +(202) 33943004 / **Mobile:** +2 (01223415321-01112221096-01011811884)

www.photoscience-co.com / info@photoscience-co.com