

Hand-Tally Differential Counter

- 5 and 8 units Lab-count with Totalizer. For simplifying differential blood cell counts.
- Represents the normal Schilling groups of Leucocytes, i.e. Baso, Eoiso., Seg., Lymph., Mono. (DBC-6) along with Myelo., Juven., Stab. (DBC-9 and DBC-8E). The above each key is the name and picture of particular cell.
- Each stroke registers both in the unit corresponding to its pressed key and in the totalizer.
- When a total of every 100 cells are counted, a bell rings or a sound buzzes. The process can be extended to 200, 300... etc.
- Counter clockwise rotate handle, it will clear all figures. Or just press "CLEAR" key to clear all figures (DBC-8E)
- DBC-8E functions CLEAR, Percentage, Decrease and Lock.

Model	DBC-6	DBC-9	DBC-8E
Leucocytes	5 units	8 units	8 units
Control	Mechanism	Mechanism	Microprocessor
Display digital	0-999	0-999	0-999
Alarm at 100	Yes	Yes	Yes
Net weight	1.1kg	1.3kg	1.5kg



DBC-8E