

Heat Stage

company offer all kinds of heat stage for the microscopes, they are with the good quality and suitable price to the customers. Not only use for the schools or the high grade lab research, these heat stages of microscopes all meet with you requests.



For Upright Biological Microscope
Room Temperature~70°C



For Inverted Biological Microscope
Room Temperature~70°C



For Polarization Microscope
Room Temperature~450°C



For Stereo Microscope
Room Temperature~70°C



Cool and Heat stage, For Biological,
Stereo Microscope
-40~50°C



For melting point test
Room Temperature~250°C



For low heat big stage
Room Temperature~250°C



For high heat big stage
Room Temperature~250°C

