


# HEATING BLOCK

## Use :

- Preservation of the samples and the reaction
- DNA tests and genetic tests increases
- Plasma (blood) clot
- The enzyme reaction, serum disabled, culture and designed to perform multiple multifunctional products.



Accurate temp. control.



Blocks by sizes



Easy Temp. control



Secured safety

## Features

- Offer the suitability with lab test by using of micro test tube, test tube, and cuvette.
- Accurate temperature controls by Digital PID controller and timer set function.
- Airtight internal design for accurate temperature control.
- Overheat protection sensor equipped. 4 different models available

## Specification

Model NO		PH-HB100	PH-HB200	PH-HB300	PH-HB400
Temp	Range	Amb 5 ~ 150 ° C			
	Controller	Digital PID Controller			
Heater		150W	240W	330W	420W
Size	Internal(WxDxH, mm)	100x80x50	100x160x50	100x235x50	200x160x50
	External(WxDxH, mm)	230x245x100	230x325x100	230x405x100	330x325x100
Safety		Over temp limiter.			
Material	Internal	Stainless steel			
	External	Steel, Powder coating			
	Plate	Aluminum			
Electric Requirements		220VAC. 50/60Hz			
Power Consumption		150W	240W	330W	420W

Order Number		Order Number	
No.	Tube Size	No.	Tube Size
1	10Ø x 20hole	5	20Ø x 12hole
2	13Ø x 20hole	6	24Ø x 6hole
3	16Ø x 12hole	7	0.5 ml x 24hole
4	18Ø x 12hole	8	1.5 ml x 24hole

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt **Code:** 12515  
**Tel. Fax:** +(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)  
[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)