



## A54.1004 Temperature Controlled Heating Stage, 20°C~50°C

A54.1004 Temperature Control Working Stage	
<b>Input Power</b>	60W , AC 110V/60HZ, or 230V/50HZ Switchable
<b>Output Power</b>	30V/50W
<b>Fuse</b>	3.15A
<b>Temp. Control Range</b>	20°C-50°C
<b>Temp. Setting Division</b>	0.1 °C
<b>Various Between Set And Display Temp.</b>	<+/-1.0°C
<b>Various Between Actual And Display Temp.</b>	<+/-1.0°C
<b>Standard Outfit</b>	Temperature Controller Box, Size 90*165*190mm - Stage: 150x130x5.8mm Cable: 1 pc Adhesive Tape: 1 pc Spare Fuse 2 pcs
<b>Packing</b>	1pc/Box, Box Size 14=17*20mm, 3 Kgs



## A54.1004 Temperature Controlled Heating Stage, 20°C~50°C



A54.1004 Temperature Control Stage Controller Box, Front View and Back View