

LOW TEMP INCUBATOR

Use : This product is used for various test such as micro organism, cultivation of tissue of plant and animal



Excellent ventilation



Parking casters



Easy to use shelves



Defrost switch



FEATURES

- Fuzzy function to prevent overheat.
- Accurate temperature controls by Digital PID controller.
- In case of temperature sensor disconnection, error message appear.
- Alarm rings when timer end.
- In case of difference between real temperature and appeared temperature, user can adjust correctly.
- When controller setting value is changed by the mistake of operation, save and recall function can restore it early.
- User can adjust temperature not to raise up certain degree from installed temperature.
- Uniform airflow by blowing fan in the chamber.
- Internal tempered glass and easy observation of the test materials by opening the external door only.
- Buffering device to enhance door adherence for insulation.
- Cable hole to use external measurement equipment.
- STS #304 used in inside chamber for durability.
- Freezer frost removal function be selected manually for the safety of samples.

SPECIFICATION

Model NO		PH-LTI81L	PH-LTI150L	PH-LTI250L
Type		1 door		
Temp	Range	-5 ° C ~ 60 ° C		
	Controller	Digital PID Controller		
Size	Internal (W x D x H, mm)	450x400x450	500x500x600	500x500x1000
	External (W x D x H, mm)	615x650x1040	665x765x1190	665x765x1590
Safety	Electrical	ELB(Earth Leakage Breaker)		
Material	Internal	Stainless steel #304		
	External	Steel, Powder coating		
	Shelves(WxD, mm)	STS #304 Wire Electro-Polished		
	Heater	Sheath heater		
	Insulation	Atiron		
	Door Packing	High Temperature grade foamed silicone rubber		
	Windows	Tempered Safety Glass		
Refrigerator		1/4HP	1/3HP	1/3HP
Electric Requirements		220VAC. 50/60Hz		
Power Consumption		900W	900W	900W

Branch:14 Fawzy Soliman St., From El Haram St., Giza, Egypt / **P.O.Box:** 127 Giza, Egypt **Code:**12515

Tel. Fax:+(202) 33943004 / **Mobile:** +2 (01223415321-01112221096-01011811884)

www.photoscience-co.com / info@photoscience-co.com