


LCD Moisture Analyzer Halogen Heating (0.001g*50g) (DSH-50-1)

- DSH-50 series,
- Moisture analyzer.
- 1mg/5mg/10mg accuracy,
- Max. range 50g.
- With RS232 interface.
- We also provide OEM service,



Model	DSH-50-1	DSH-50-5	DSH-50-10
Capacity (g)	50 g		
Readability (g)	0.001 g	0.005 g	0.01 g
Repeatability (g)	0.2 %	0.3 %	0.5 %
Min. Weight of Sample (g)	0.5 g		
Propose Sample Weight (g)	0-10 g		
Heating Time (M.)	1 – 99 minutes (stepping is 1 m.)		
Heating	Standard		
Working Temperatuer (Centi Degree)	10 - 30 (Centi Degree)		
Interface	Rs232		
Time Control	Setting time by artificially or automatically		
Heating Temperatuer (Centi Degree)	50 - 160 Centi Degree (stepping is 1)		
Display	Moisture %, Solid %, Weight, Time, Temperature		
Pan (mm)	100 mm		
Outline (mm)	285 × 160 × 150 mm		
Shipping Size (mm)	455*328*350 mm		
N. W. / G. W. (kg)	3 kg / 4.2 kg		
Heating Way	Halogen		

**Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt P.O.Box: 127 Giza, Egypt Code:12515
Tel. Fax:+(202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)**

www.photoscience-co.com / info@photoscience-co.com