

# CIRCUCLATOR OIL BATH

**Use :** Circulator Oil Bath is generally used for basic laboratory equipments for biotech and industrial fields. For test reliability and ease of use, this products has following characteristics;



## Specification

Model NO		PH-COBD7L	PH-COBD12L	PH-COBD25L
Type		Circulator Oil Bath		
Controller		Digital PID Controller		
Range		Max 200□		
Display		LED 4 Digital Display		
Temp	Accuracy	±0.2□		
Size	Internal (WxDxH, mm)	300x160x150	330x220x170	400x350x180
Circulator Pump		Circulator Pump		
Safety		Over temp limiter.		
Material	Internal	STS #304, Stainless steel		
	External	Steel, Powder coating		
Electric Requirements		220VAC 50Hz or 60Hz		
Heater		1500W	2000W	3000W

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt **Code:**12515  
**Tel. Fax:**+(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)  
[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)