

# SHKER

Forward / Reverse Configurable Orbital Shaker.  
Slim and can be used on a table in a compact design.  
Minimized vibration of the product to a low center of gravity design,  
Durability is excellent in stir continued use of the belt and flexible.

## Feature

And durability, this ability is excellent loading capacity  
A low center of gravity to minimize vibration in the product design.  
It said flexible belt is applied to prevent factors were often broken.  
Direction of the shaker is possible to forward rotation and reverse rotation and in accordance with the time set by adding a timer function  
Shaker operation will be written.  
Depending on the weight of the sample was designed, which is not affected by the gravity center position.  
RPM changes are possible during operation, and not constrained in voltage.  
Sticky meth alcohol, it is possible to use a semi-permanent and washed with water



## SPECIFICATION

MODEL	PURIKER 350
Dimension_IN(WxDxH)mm	350 x 395 x 93
Dimension_Aluminium platform capacity	372 x 372 x 15
Dimension_Aluminium spring platform capacity	375 x 375 x 75 spring 64EA
100ML	25
250ML	16
500ML	9
1000ML	5
2000ML	4
Material (exterior)	Steel plate with powder coating
Controller	Microprocessor control
Speed.Range	30 ~ 350 rpm
Speed.increment	10 rpm
Time	Range Selectable 59m.59s, 59h.59m, continuous
Motor	Brushless DC motor
Motor driving direction	Forward, reverse, auto reverse
Stroke	22mm
Display	7-Segment LED display
Operating input	Touch botton
Power	European 220 VAC 50/60Hz STD 1 phase

**Branch:** 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515

**Tel. Fax:** +(202) 33943004 / **Mobile:** +2 / 01223415321-0111221096-01011811884

[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)