

SHAKING INCUBATOR HORIZONTAL

Use : This product is used as the general and basic laboratory equipment such as industrial and Biotechnology field



Parking Caster



Multi Flask holder plate



Clifton lamp



Membrane touch pad controller



Features

- Accurate temperature controls by Digital PID controller.
- Accurate speed control and digital mark function made for user convenience.
- Circulation device in internal chamber have a little noise and keep overall temperature evenly.
- Overheat, over current prevention device for users and safety function to stop shaking device automatically when door open.
- Various holder and plate available according to users spec.
- Little noise designed structurally

Specification

Model NO		PH-SIH300R	PH-SIH300RF
Temp	Range	0 °C ~ 60 °C	
	Controller	Digital PID Controller	
	Accuracy	±0.2 °C at 37 °C	
Size	Internal (W x D x H, mm)	890 x 640 x 400	
	External (W x D x H, mm)	1120 x 720 x 940	
Safety	Electrical	Over temp limiter, ELB(Earth Leakage Breaker), Pause stop	
Material	Internal	Stainless steel STS #304	
	External	Steel, Powder coating	
Speed Range		20 ~300 rpm	
Shaking Stroke		20mm Orbital Motion	
Light System		None	FL 20W * 4ea
Refrigerator		1/3 HP	1/3 HP
Electric Requirements		220VAC. 50/60Hz	

Clamp Numbers	Plate Size (mm)
Plate Size(mm)	810×560
100 ml	77ea
250 ml	48ea
500 ml	24~25ea
1000 ml	20ea
2000 ml	12ea
3000 ml	11ea

Max load (Kg)	Model	PH-SIH300R	PH-SIH300RF
		Kg	80

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt / **P.O.Box:** 127 Giza, Egypt **Code:**12515
Tel. Fax: +(202) 33943004 / **Mobile:** +2 (01223415321-01112221096-01011811884)

www.photoscience-co.com / info@photoscience-co.com