

SHAKING WATER BATH DIGITAL

Use : This product is used as the general and basic laboratory equipment such as industrial and Biotechnology field.



Stable hanger



Low noise and multi rack



Attachable multi rack



Features

- The grounded circuit prevent a short circuit and over current of electricity.
- Time set function and constant temperature.
- Unique design multi rack separable (utility model)
- Minimize noise to solve the problem of old products.
- Equipped internal bottom for the safety of the user and protection of the unit.

Specification

Model NO		PH-SWBD25L	PH-SWBD35L	PH-SWBD55L
Type		Digital		
Controller		Digital PID Controller		
Range		Ambient +5~100°C		
Temp	Accuracy	±0.2°C		
Size	Internal(W x D x H, mm)	350x300x240	500x300x240	660x350x240
	External(W x D x H, mm)	580x360x400	730x360x400	890x410x400
Safety		Over Temp limiter.		
Material	Internal	Stainless steel STS #304		
	External	Steel, Powder coating		
External		220VAC, 50 or 60Hz		
Rack		Spring wire rack		
Heater		1500W	2300W	2800W

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt **Code:**12515
Tel. Fax: +(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)
www.photoscience-co.com / info@photoscience-co.com