

VISIBLE BATH DIGITAL


Use : This product is used as the general and basic laboratory equipment such as industrial and Biotechnology field.



Embedded handle



Viscosity install hall



Circulation hose



Tempered glass



Features

- Equipped overheat protection sensor.
- Time set function and constant temperature.
- Uniquely designed smooth exterior and mono piece bath.
- Equipped internal bottom for the safety of the user and protection of the unit.

Specification

Model NO	PH-VBD30L	PH-VBD55L
Camber volume	30l	55l
Range	Amb +5°C~100°C	
Controller	Digital PID Controller	
Internal(W x D x H, mm)	370x220x370	500x300x370
External(W x D x H, mm)	450x285x715	570x365x715
Accuracy	±0.2°C	
Safety	Over temp limiter.	
Internal	Stainless steel #304	
External	Steel, Powder coating	
Heater	Sheath Heater	
Windows	Tempered Safety Glass	
Electric Requirements	220VAC. 50/60Hz	
Power Consumption	2400W	4000W

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515

Tel. Fax: +(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)

www.photoscience-co.com / info@photoscience-co.com