

CIRCULATOR WATER BATH

Use : This product is used as the general and basic laboratory equipment such as industrial and Biotechnology field.




Circulation hose



PID controller



Drainage line



Water bath cover



Features

- The grounded circuit prevent a short circuit and over current of electricity.
- Anti-corrosive by STS #304 and efficient insulation structure.
- External circulation function by strong circulation pump
- Time set function and constant temperature.
- Uniquely designed smooth exterior and circulation pump separable

Specification

Model NO		PH-CWBD13L	PH-CWBD20L	PH-CWBD30L
Type		Digital		
Controller		Digital PID Controller		
Range		Amb +5°C~100°C		
Temp	Accuracy	±0.2□		
Size	Internal(WxDxH, mm)	340x210x180	400x290x200	450x300x230
Safety		Over temp limiter.		
Material	Internal	Stainless steel STS #304		
	External	Steel, Powder coating		
Electric Requirements		220VAC 50/60Hz		
Heater		1500W	2000W	2000W
Cover		Stainless steel #304		

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt P.O.Box: 127 Giza, Egypt Code:12515
Tel. Fax:+(202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)
www.photoscience-co.com / info@photoscience-co.com