


WATER BATH

Use : This product is used as the general and basic laboratory equipment such as industrial and Biotechnology field.



Features

- Equipped overheat protection sensor.
- Time set function and constant temperature.
- Uniquely designed smooth exterior and mono piece bath.
- Equipped internal bottom for the safety of the user and protection of the unit.

Specification

Model NO		PH-WBD06	PH-WBD11	PH-WBD23
Type		Digital		
Controller		Digital PID Controller		
Range		Ambient +5~100° C		
Temp	Accuracy	±0.2° C		
Size	Internal (WxDxH, mm)	240x140x150	300x240x150	515x330x150
	External (WxDxH, mm)	320x270x220	360x370x230	580x410x230
Safety		Over temp limiter.		
Material	Internal	Seamless Stainless steel		
	External	Steel, Powder coating		
Electric Requirements		220VAC 50/60Hz		
Heater		500W	1000W	2400W

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt **Code:**12515
Tel. Fax:+(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)
www.photoscience-co.com / info@photoscience-co.com