

El destilador de vidrio BASIC PH4 está especialmente diseñado para la destilación en continuo de un modo muy sencillo, con total seguridad y confianza gracias al termostato de seguridad incorporado de que dispone.

Cuenta con una resistencia de cobre recubierta de níquel-cromo de una potencia de 3Kw que permite una producción de agua destilada de 4 litros a la hora, incluida en un depósito de destilación fabricado en vidrio borosilicato de alta calidad.

El equipo se suministra con marcado y certificado CE.

The glass water still BASIC PH4 is specially design for an easy and continuous distillation process with total safety and confidence thanks to its built-in safety thermostat.

Provided with a 3 Kw nickel-chromium plated copper resistance, the water still BASIC PH4 has a production of 4 L/h of distilled water.

The horizontal boiler is made of high quality borosilicate glass. CE certified.



Conductividad (Agua corriente) <i>Conductivity (with tap water)</i>	3.0-4.0 $\mu$ S/cm
Conductividad (Agua tratada) <i>Conductivity (with treated water)</i>	1.5-2.0 $\mu$ S/cm
pH	5.5-6.0
Dimensiones / Dimensions	50x15x45 cm

## Accesorios / Accessories



En caso de deterioro, roturas accidentales etc, de las partes que constituyen el equipo, están disponibles a continuación los siguientes componentes para su reposición

In case of damage, breakage, etc. of any of the components of the equipment; there are available the following pieces for replacement.



Referencia Code	Descripción Description
99002100	Bloque calefactor, con resistencia / Heating block, with resistance
99002104	Depósito de vidrio / Glass boiler
99002105	Refrigerante de vidrio / Glass condenser
99002103	Juego de mangueras con conectores rojos / Hose kit with red connectors
99002101	Juntas+rosca+2 sujeciones al destilador / Joint hit+screw+2 fixers for boiler
99002107	Soporte metálico / Metal support
99002108	2 sujeciones al destilador / 2 fixers for boiler